



Medical Grade Silicone Biliary Stents

Our Biliary Stents are manufactured from high quality, medical grade silicone material. The soft consistency of our stents provides for patient comfort while the biocompatibility of our medical grade silicone material supports short term placement. These stents are easily placed percutaneously and can be changed readily.

As shown in the illustration below, each stent contains 33 side-holes for efficient drainage. A radiopaque tantalum marker is placed in one of the most proximal side-holes to assist confirmation of correct placement after insertion. Our Silicone Biliary Stents are available in a variety of sizes from 8F to 29F; all are 24 inches in length but can be cut to adapt to the anatomy of the individual patient.



Medical grade plastic luer lock adaptors are included with the 8F to 22F stents to facilitate ease of connection to an external drainage system, if needed. The Bentec Medical Silicone Biliary Stents ship sterile.

Medical Grade Silicone Biliary Stents

FEATURES	BENEFITS
<ul style="list-style-type: none"> Medical grade silicone – soft, flexible design Radiopaque tantalum marker in proximal side-hole 33 side-holes on each stent Range of sizes from 8F to 29F Luer lock adaptor included with stents sizes 8F to 22F Sterile product Not made with natural rubber latex 	<ul style="list-style-type: none"> Biocompatibility supports short term placement (29 days or less) Decreases potential patient discomfort Enables confirmation of acceptable stent location by x-ray post insertion Provides efficient drainage capability Provides options to address different patient anatomies & drainage requirements Facilitates connection to percutaneous drainage systems Reduces potential of infection Eliminates potential for latex-induced allergic reactions



Ordering Information

SILICONE BILIARY STENTS						STERILE
Product Code	Size (F)	Length (inches)	Radiopaque Marker	Luer Adaptor	Quantity	
GR83064-12	8	24	YES	YES	1 per pkg	
GR83064-03	10	24	YES	YES	1 per pkg	
GR83064-01	12	24	YES	YES	1 per pkg	
GR83064-02	14	24	YES	YES	1 per pkg	
GR83064-04	16	24	YES	YES	1 per pkg	
GR83064-05	18	24	YES	YES	1 per pkg	
GR83064-08	19	24	YES	YES	1 per pkg	
GR83064-06	20	24	YES	YES	1 per pkg	
GR83064-07	22	24	YES	YES	1 per pkg	
GR83064-10	24	24	YES	NO	1 per pkg	
GR83064-13	26	24	YES	NO	1 per pkg	
GR83064-14	28	24	YES	NO	1 per pkg	
GR83064-11	29	24	YES	NO	1 per pkg	

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
 All Bentec Medical products are manufactured in the U.S. by Bentec Medical in compliance with strict CGMP Guidelines.
 Copyright 2017 Bentec Medical.

PDS0001.A July 2017