



Silicone Bulk Tubing

Bentec Medical's platinum-cured Silicone Bulk Tubing is translucent, smooth, soft and pliable. It is available in a range of sizes as shown in the table below. Silicone Bulk Tubing is available in 50 foot lengths only and ships non-sterile.

Ordering Information

SILICONE BULK TUBING				NON-STERILE
Product Code	ID (in)	OD (in)	Length (ft)	Quantity
BM3-1050	.012	.025	50	1 per pkg
BM3-1350	.020	.037	50	1 per pkg
BM3-1550	.025	.047	50	1 per pkg
BM3-1750	.030	.065	50	1 per pkg
BM3-1850	.031	.156	50	1 per pkg
BM3-2050	.040	.085	50	1 per pkg
BM3-2350	.058	.077	50	1 per pkg
BM3-2650	.062	.095	50	1 per pkg
BM3-2850	.062	.125	50	1 per pkg
BM3-2950	.062	.188	50	1 per pkg
BM3-3050	.078	.125	50	1 per pkg
BM3-3250	.104	.192	50	1 per pkg
BM3-3350	.132	.183	50	1 per pkg
BM3-3650	.125	.250	50	1 per pkg
BM3-4050	.187	.312	50	1 per pkg
BM3-4250	.187	.375	50	1 per pkg
BM3-4450	.250	.375	50	1 per pkg