


Medical Grade Silicone Blocks

Bentec Silicone Blocks are rectangular-shaped pieces of sterile implant-grade elastomeric silicone. The dimensional size of our product varies in length, width, and thickness as well as the consistency of the Silicone. Bentec Silicone Blocks are available in four different consistencies: Soft, Medium, Firm and Extra-Firm with durometers ranging from Shore 20A to 65A.

The Blocks may be readily trimmed with a knife or scissors to the desired configuration. This feature allows surgeons to custom fabricate an implantable prosthesis for their specific indication during surgery.



Medical Grade Silicone Blocks

FEATURES	BENEFITS
<ul style="list-style-type: none"> Implant-grade silicone elastomeric Trimmable with knife or scissors Range of block size & silicone consistency Sterile product Not made with natural rubber latex 	<ul style="list-style-type: none"> Ensures product quality & performance for either long term or short term implantation Provides surgeon with flexibility to adapt product during surgery to specific anatomical or procedural need at time of surgery Provides surgeon with options to address individual procedural needs Reduces potential of infection Eliminates potential for latex-induced allergic reactions



Ordering Information

SILICONE BLOCKS			STERILE
Product Code	Size	Consistency	Quantity
PR88032-222	1 in x 1 in x 0.5 in	Medium	1 per pkg
PR88032-31	1.91cm x 1.91cm x .95cm	Soft	1 per pkg
PR88032-32	1.91cm x 1.91cm x .95cm	Medium	1 per pkg
PR88032-33	1.91cm x 1.91cm x .95cm	Firm	1 per pkg
PR88032-34	1.91cm x 1.91cm x .95cm	X-Firm	1 per pkg
PR88032-21	7cm x 4cm x 1.5cm	Soft	1 per pkg
PR88032-22	7cm x 4cm x 1.5cm	Medium	1 per pkg
PR88032-23	7cm x 4cm x 1.5cm	Firm	1 per pkg
PR88032-24	7cm x 4cm x 1.5cm	X-Firm	1 per pkg
PR88032-01	7cm x 12cm x 1.5cm	Soft	1 per pkg
PR88032-02	7cm x 12cm x 1.5cm	Medium	1 per pkg
PR88032-03	7cm x 12cm x 1.5cm	Firm	1 per pkg
PR88032-04	7cm x 12cm x 1.5cm	X-Firm	1 per pkg

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
 All Bentec Medical products are manufactured in the U.S. by Bentec Medical in compliance with strict CGMP Guidelines.
 Copyright 2017 Bentec Medical.

PDS0004.B September 2017