



Medical Grade Silicone Sheeting

Our sterile, implantable grade, Silicone Sheeting is designed for either long or short term applications in or on the body. Silicone Sheets are available from Bentec Medical in a wide range of sizes & thicknesses. These sheets can be shaped by the surgeon into different configurations & sizes to meet individual procedure requirements providing flexibility to customize for each patient.

Bentec Medical Silicone Sheeting is selected by surgeons for many different procedures, including the construction of "chimneys" for neonates with gastroschisis or omphaloceles*, reinforcement of wound or surgical incision closures, scar reduction and reconstructive plastic surgery such as abdominoplasty, breast reduction and breast augmentation.

*Ask your Bentec Medical representative about the Bentec Medical Spring-Loaded Silicone Silo Bag (the 'Gold Standard' for staged reductions of ventral wall defects).

FEATURES	BENEFITS
<ul style="list-style-type: none"> ▪ Mesh Reinforced, Non-reinforced & Open Border options 	<ul style="list-style-type: none"> ▪ Mesh Reinforced sheets minimize stretching while allowing suture placement ▪ Provides surgeon with multiple options to address individual procedural needs
<ul style="list-style-type: none"> ▪ Medical grade silicone 	<ul style="list-style-type: none"> ▪ Quality in both long term & short term applications
<ul style="list-style-type: none"> ▪ Translucence 	<ul style="list-style-type: none"> ▪ Allows visualization of tissue post application
<ul style="list-style-type: none"> ▪ Range of thickness: 0.005" to 0.125" ▪ Standard Sizes (3"x 4" & 6"x 8") with additional sizes up to 14"x 18" available 	<ul style="list-style-type: none"> ▪ Provides surgeon with selection of sheeting thickness & size to address individual procedural needs
<ul style="list-style-type: none"> ▪ Sterile product 	<ul style="list-style-type: none"> ▪ Reduces potential of infection
<ul style="list-style-type: none"> ▪ Not made with natural rubber latex 	<ul style="list-style-type: none"> ▪ Eliminates potential for latex-induced allergic reactions

Silicone Sheeting

Ordering Information

Silicone Sheeting				STERILE			
<i>Non-Reinforced</i>				<i>Reinforced</i>			
Product Code	Size (inches)	Thickness (inches)	Quantity	Product Code	Size (inches)	Thickness (inches)	Quantity
PR72034-66N	1 X 1	0.005	2 per pkg	PR72034-51R	2 X 3	0.010	2 per pkg
PR72034-62N	2 X 2	0.020	2 per pkg	PR72034-52R	2 X 3	0.020	2 per pkg
PR72034-64N	2 X 2	0.040	2 per pkg	PR72034-53R	2 X 3	0.030	2 per pkg
PR72034-51N	2 X 3	0.010	2 per pkg	PR72034-54R	2 X 3	0.040	2 per pkg
PR72034-52N	2 X 3	0.020	2 per pkg	PR72034-55R	2X3	0.050	2 per pkg
PR72034-53N	2 X 3	0.030	2 per pkg	PR72034-007R	6 x 8	0.007	1 per pkg
PR72034-54N	2 X 3	0.040	2 per pkg	PR72034-01R	6 x 8	0.010	1 per pkg
PR72034-20N	3 x 4	0.005	2 per pkg	PR72034-02R	6 x 8	0.020	1 per pkg
PR72034-21N	3 x 4	0.010	2 per pkg	PR72034-03R	6 x 8	0.030	1 per pkg
PR72034-22N	3 x 4	0.020	2 per pkg	PR72034-04R	6 x 8	0.040	1 per pkg
PR72034-23N	3 x 4	0.030	2 per pkg	PR72034-05R	6 x 8	0.050	1 per pkg
PR72034-24N	3 x 4	0.040	2 per pkg	PR72034-06R	6 x 8	0.060	1 per pkg
PR72034-00N	6 x 8	0.005	1 per pkg	PR72034-07R	6 x 8	0.070	1 per pkg
PR72034-01N	6 x 8	0.010	1 per pkg	PR72034-08R	6 x 8	0.080	1 per pkg
PR72034-02N	6 x 8	0.020	1 per pkg	PR72034-10R	6 x 8	0.100	1 per pkg
PR72034-03N	6 x 8	0.030	1 per pkg	PR72034-12R	6 x 8	0.125	1 per pkg
PR72034-04N	6 x 8	0.040	1 per pkg	PR72034-60R	2 x 2	0.007	2 per pkg
PR72034-05N	6 x 8	0.050	1 per pkg	PR72034-63R	2 x 2	0.030	2 per pkg
PR72034-06N	6 x 8	0.060	1 per pkg	PR72034-64R	2 x 2	0.040	2 per pkg
PR72034-07N	6 x 8	0.070	1 per pkg	PR72034-62R	2 x 2	0.020	2 per pkg
PR72034-08N	6 x 8	0.080	1 per pkg	PR72034-14R	6 x 8	0.140	1 per pkg
PR72034-09N	6 x 8	0.090	1 per pkg	PR72034-407R	14 x 18	0.007	1 per pkg
PR72034-10N	6 x 8	0.100	1 per pkg				
PR72034-11N	6 x 8	0.120	1 per pkg				
PR72034-12N	6 x 8	0.125	1 per pkg				



Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
 All Bentec Medical products are manufactured in the U.S. by Bentec Medical in compliance with strict GMP Guidelines.
 Copyright 2017 Bentec Medical

PDS0005.A September 2017