


Medical Grade Silicone Elastomer Tubing



Bentec Medical's Silicone Elastomer Tubing is manufactured from the highest quality medical grade Silastic™ silicone. It is smooth, soft and pliable, making it easy to use and comfortable for the patient. Our Silicone Elastomer Tubing may be used in a variety of surgical applications & procedures based on clinician preference. Bentec Medical's sterile Silicone Elastomer Tubing is available in 2 configurations (Translucent or Radiopaque), a wide range of French sizes and a standard length of 100 cm (39.37 inches).

FEATURES

- Smooth, soft and pliable material
- 2 Configurations – Translucent or Radiopaque
- Range of tubing sizes
- Sterile product
- Not made with natural rubber latex

BENEFITS

- Easy to use for clinician
- Provides patient comfort when in place
- Translucent enables visualization of contents in the tubing
- Radiopaque enables visualization of placement/location in the body via X-ray
- Usable for a variety of applications in different procedures based on clinician need or preference
- Reduces potential for infection
- Eliminates potential for latex-induced allergic reactions

Silicone Elastomer Tubing

Ordering Information

Silicone Tubing						STERILE
Product Codes		Size (Fr)	OD (in)	ID (in)	Length (in)	Quantity
Translucent	Radiopaque*					
GS75160-02	GS75160-42	2	.025	.012	39.37	1 per pkg
GS75160-03	GS75160-43	3	.037	.020	39.37	1 per pkg
GS75160-04	GS75160-44	4	.052	.027	39.37	1 per pkg
GS75160-05	GS75160-45	5	.065	.030	39.37	1 per pkg
GS75160-06	GS75160-46	6	.078	.047	39.37	1 per pkg
GS75160-07	GS75160-48	7	.092	.055	39.37	1 per pkg
GS75160-08	GS75160-49	8	.105	.063	39.37	1 per pkg
GS75160-09	GS75160-50	9	.118	.071	39.37	1 per pkg
GS75160-10	GS75160-51	10	.131	.078	39.37	1 per pkg
GS75160-11	GS75160-52	11	.142	.110	39.37	1 per pkg
GS75160-12	GS75160-54	12	.155	.105	39.37	1 per pkg
GS75160-14	GS75160-55	14	.183	.122	39.37	1 per pkg
GS75160-15	GS75160-56	15	.197	.133	39.37	1 per pkg
GS75160-16	GS75160-58	16	.210	.140	39.37	1 per pkg
GS75160-18	GS75160-59	18	.240	.160	39.37	1 per pkg
GS75160-19	GS75160-60	19	.250	.133	39.37	1 per pkg
GS75160-20	GS75160-62	20	.262	.174	39.37	1 per pkg
GS75160-22	GS75160-64	22	.290	.200	39.37	1 per pkg
GS75160-24		24	.314	.209	39.37	1 per pkg
GS75160-26		26	.335	.295	39.37	1 per pkg
GS75160-28		28	.360	.181	39.37	1 per pkg

*Radiopaque Silicone Elastomer Tubing in 4, 8 and 12 inch lengths available by special order. Ask your Bentec Medical representative for more information.

