



Medical Grade Silicone Tracheostoma Vents


Our Tracheostoma Vent is designed to maintain the tracheostomy stoma following tracheal and/or laryngeal procedures when a full size tracheostomy tube is no longer necessary. The Bentec Medical Tracheostoma Vent is constructed of soft, flexible, medical grade silicone and our unique design configuration maintains placement securely without the need for neck strap or ties. Available in various sizes & configurations, our Tracheostoma Vent provides a patient friendly option for the long-term maintenance of a tracheostomy stoma.

FEATURES	BENEFITS
<ul style="list-style-type: none"> ▪ Medical grade silicone – soft, flexible material ▪ Unique design ▪ Range of sizes & configurations ▪ Not made with natural rubber latex 	<ul style="list-style-type: none"> ▪ Helps to <ul style="list-style-type: none"> – reduce problems with stoma opening scarring – provide more comfort in extreme neck positions – reduce risk of tracheal wall perforations ▪ Biocompatible ▪ Provides for a snug, secure fit without a neck strap ▪ Supports easy removal & cleaning if needed ▪ Fits a variety of stoma sizes & patient anatomies ▪ Allows enlargement of stoma breathing area over time if needed ▪ Eliminates potential for latex induced allergic reactions

Tracheostoma Vents

Configuration

Standard Vents



A = OD in mm B = ID in mm C = Length in mm

Ordering Information

Tracheostoma Vent

NON-STERILE

Product Code	Size	OD (mm)	ID (mm)	Length (mm)	Quantity
530-7005	5 standard	9.0	6.2	14.5	1 per pkg
530-7006	6 standard	9.5	8.0	15.3	1 per pkg
530-7007	7 standard	11.0	9.5	16.5	1 per pkg
530-7008	8 standard	12.0	9.5	18.0	1 per pkg
530-7009	9 standard	13.5	10.5	19.0	1 per pkg
530-7010	10 standard	15.0	12.0	20.5	1 per pkg
530-7012	12 standard	17.0	13.5	20.5	1 per pkg



Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
All Bentec Medical products are manufactured in the U.S. by Bentec Medical in compliance with strict GMP Guidelines.
Copyright 2017 Bentec Medical.

PDS0009.B September 2017